

	<b>Test Protocol</b>	
--	----------------------	--

Date: \_\_\_\_\_  
Testing Personnel: \_\_\_\_\_

Automatic Storage System		Biobank Management System	
System		System	
Software version		Software version	
Consumables		Consumables	

No.	Test Case	na	Passed?	AAS	BIMS	Comments
1	All samples available and in the right place and order	<input type="checkbox"/>				
2	Sample / Rack ID early in order	<input type="checkbox"/>				
3	Sample / Rack ID late in order	<input type="checkbox"/>				
4	Sample/ rack is missing	<input type="checkbox"/>				
5	Too many samples/ racks	<input type="checkbox"/>				
6	Sample / Rack ID unknown	<input type="checkbox"/>				
7	Samples/ Racks of the orders are swapped	<input type="checkbox"/>				
8	No order available	<input type="checkbox"/>				
9	Sample/ rack size too small	<input type="checkbox"/>				
10	Sample/ rack size too big	<input type="checkbox"/>				